Distributed by :

Peters SURGICAL

42 RUE BENOIT FRACHON - 93303 BOBIGNY - France
TEL: 33 (0) 1 48 10 62 62 - FAX : 33 (0) 1 48 91 22 99
E-MAIL : peters@peters-surgical.com
www.peters-surgical.com

PROMEPLA S.A.M

9 avenue Albert II - 98000 MONACO

REF

CAN2505

CAN25010

CAN2505GL

CAN2505P

CAN2505DP

CAN2505W

€ 0459

Fr

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INSTRUCTIONS FOR USE SINGLE-USE STERILE IRRIGATIONSUCTION KIT

En

CAUTION

Please read these Instructions carefully before using the single-use irrigation-suction kit with cannula. Failure to comply with the following recommendations may entail serious consequences. Laparoscopic surgery should only be performed by qualified doctors having been specifically trained for the purpose and having gained sufficient experience of use. All surgeons performing laparoscopic operations should be aware of the possible risks and complications involved using such surgical technique (notably risks for the patient and in manipulating instruments connected to an external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure). Ensure that the patient has been duly informed of the operation procedure, such as required by legislation in force.

DESCRIPTION OF THE INSTRUMENT

The single-use irrigation-suction kit with cannula, comprised of a double piston handpiece and a polytetrafluoroethylene (PTFE)-covered cannula, is proposed in the following variants:

References	Cannula Diameter	Cannula Length	Tube Diameter	Tube Length	Percussor	Distinctive Features
CAN2505	5 mm	330 mm	5 mm	2 500 mm	1	Kit with standard-length cannula.
CAN2505GL	5 mm	450 mm	5 mm	2 500 mm	1	Kit with extra-long cannula.
CAN2505W	5 mm	330 mm	5 mm	2 500 mm	1	Kit with standard-length cannula and pump adapter tip (8 mm internal diameter on pump end; 4 mm internal diameter on percussor end).
CAN2505P	5 mm	330 mm	5 mm	2 500 mm	1	Kit with standard-length cannula, anti-reflux system.
CAN2505DP	5 mm	330 mm	5 mm	2 500 mm	2	-
CAN25010	10 mm	330 mm	5 mm	2 500 mm	1	Kit with large-diameter cannula.

A cap is fitted at the rear of the double-piston handpiece.

The irrigation line can be seen on the blue piston, with its tube being equipped with a clamp and a pouch percursor (except CAN2505DP)

The irrigation line can be seen on the red piston, with its tube being equipped with a flexible universal connector

↳ **On no account may Promepla be held liable in the case of failure by the operator to observe the recommendations pertaining to proper surgical practices.**

INDICATIONS OF USE

The irrigation-suction kit with cannula is a single-use device for laparoscopic surgery, serving to cleanse the abdominal cavity and to suck up any waste. It should only be used by qualified practitioners having gained due experience.

CONTRAINDICATIONS

Use of the single-use irrigation-suction kit with cannula is strictly forbidden in the case of any contraindication to undergo laparoscopic surgery. On no account may Promepla be held liable for the assessment of the contraindication in relation to laparoscopic surgery.

WARNINGS

↳ **Specific warnings with regard to surgery are stated herein, throughout the present Instructions (in bold type and preceded by L).**

↳ **These Instructions solely contain information on the irrigation-suction kit, as well as on how to use it. On no account should they be considered as an instructions' manual for using irrigation-suction kits in general.**

↳ **Irrigation-suction kits are designed for single use only. On no account may they be used a second time over, whatever the reason or the means.**

↳ **The irrigation-suction kits are delivered in sterile packaging. They should never be resterilized, regardless of whether they have been used or not.**

↳ **Any reuse or resterilization of this instrument shall affect its original sterility and may lead to infection or to cross-contamination in patients. Reuse or resterilization may also compromise the integrity and functioning of the instrument, possibly causing injury, sickness or even death of a patient.**

↳ **Responsibility with regard to the storage, the opening of the packaging, the appropriate use, handling, insertion and the destruction of the irrigation-suction kits lies solely with the operator. On no account may Promepla be held liable in the case of failure by the operator to observe proper surgical practices for the storage, the opening, the appropriate use, the manipulation, the insertion and the destruction of the irrigation-suction kits. Likewise, the manufacturer may not be held liable in the case of reuse and/or resterilization.**

↳ **The sterile packaging, along with its contents, should be destroyed in any of the following cases:**

- Packaging severely crushed or clearly damaged;
- Visible damp marks on the cardboard packaging;
- Perforation of the plastic packet.

STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment adapted to sterile medical devices.
- For transportation, the above points should be observed.

Information for the United States (USA):

- US Federal Law only allows the sale of this device by a doctor or upon presentation of a doctor's prescription.
- Sterilized with ethylene oxide.
- Single use only - Not to be reused.

INSTRUCTIONS FOR USE

Before insertion of the single-use irrigation-suction kit, the operation zone, i.e. the abdominal wall, should be prepared as usual, according to proper surgical practices.

Prior to using the single-use irrigation-suction kit with cannula, please ensure that such cannula is correctly adapted to the diameters of the trocars being used.

1. Open the multi-pack box and take out the relevant medical device.
2. Make sure that the sterile packaging protecting the medical device is neither torn nor damaged. If in any doubt, however slight, destroy the sterile packaging, together with its contents.
3. Check the use-by date printed on the packaging.
4. Open the sterile peel-back packaging in accordance with the proper practices for operation zones.

↳ **Irrigation-suction kits and their packaging should be destroyed as of when sterility may no longer be guaranteed, notably in the following cases: expiry of the use-by date on the packaging; packaging severely crushed or clearly damaged; damp marks visible on the packaging (both the cardboard box and the plastic packet); perforation of the plastic packet.**

5. Check:
 - the good state of the irrigation-suction kit;
 - the existence of all components: tubes, cannula, handpiece cap;
 - that the cannula and the rear cap are tightly screwed.

↳ **Damaged irrigation-suction kits should not be used.**

6. Remove the cap from the cannula (except CAN2505P)
7. Retain the irrigation cannula, the handpiece and approximately one metre of tube within the operation zone.
8. Before starting surgery, make sure that all the laparoscopic instruments are adapted to the single-use irrigation-suction kits.

↳ **On no account may Promepla be held liable in the case of failure by the operator to observe proper surgical practices when verifying whether the laparoscopic instruments are compatible with irrigation-suction kits.**

↳ **Be sure to follow all regular techniques according to proper surgical practices in order to protect the vessels and organs from any possible injury.**

↳ **CAUTION: For CAN2505P, reduction in the internal section of the cannula may reduce the suction of any large clots (or any significant waste).**

9. Shut the clamp of the irrigation tube and, after removing the percursor's protective cap, percuss the pouch or the flask of solution in a sterile manner. Where necessary, use an adaptor for connection to a pump (only CAN2505W).
10. Connect the flexible connector of the suction tube to the main vacuum unit via retention jars. The body liquids and/or cleansing solutions are sucked out via the main vacuum unit within the operating theatre.

11. Prior to using the cannula, open the clamp of the irrigation tube in order to allow the liquid to pass through and to flush out any air inside the tube.

12. During the course of use:
 - the bleu piston serves for irrigation, while the red piston serves for suction;
 - it is recommended to only release the red suction piston once all the liquid contained in the cannula and the handpiece has been sucked out;

- avoid jolting the piston whenever activated.

- During a long operation or when the aspirated fluids are particularly polluted, it is recommended to perform a flush by pressing both pistons simultaneously for several seconds.

In order to change the pouch, shut the irrigation tube clamp and replace the pouch or the flask of solution in a sterile manner.

At the end of surgery, shut the irrigation tube clamp, disconnect the flexible universal suction connector and duly dispose of the cannula-handpiece-tube kit.

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