
# SonoGuard

PROMECON

- DE** Sterile Abdeckung für Ultraschallsonden
- EN** Sterile cover for ultrasound probes
- FR** Housse stérile pour sondes à ultrason
- IT** Copertura sterile per sonda ad ultrasuoni
- ES** Funda estéril para sondas de ecografía
- PT** Capa de proteção esterilizada para sondas de ultrassom
- NL** Steriele hoes voor ultrasone probes
- DK** Steril afdækning til ultralydsprober
- FI** Steriili suojus ultraäänisondeihin
- NO** Steril tildekning til ultralydssonder
- SE** Sterilt skydd för ultraljudssonder

- |  |                             |
|--|-----------------------------|
| <b>DE</b> GEBRAUCHSANWEISUNG           | <b>NL</b> GEBUIKSAANWIJZING |
| <b>EN</b> INSTRUCTIONS FOR USE         | <b>DK</b> BRUGSANVISNING    |
| <b>FR</b> MODE D'EMPLOI                | <b>FI</b> KÄYTTÖOHJE        |
| <b>IT</b> ISTRUZIONI PER L'USO         | <b>NO</b> BRUKSANVISNING    |
| <b>ES</b> INSTRUCCIONES DE UTILIZACIÓN | <b>SE</b> ANVÄNDARANVISNING |
| <b>PT</b> INSTRUÇÕES DE UTILIZAÇÃO     |                             |

**CE**  
**0537**

	<p>Promecon GmbH          Kreuzlerstrasse 10          20095 Hamburg          Germany          Tel.: +49 40 3690 1690          Email: info@promecon-medical.de          www.promecon-medical.com</p>
<p><b>STERILE EO</b> <b>LATEX-FREE</b> <b>DEHP-FREE</b></p>	

### **Wichtig**

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**Important**

Please read the following information and instructions for use carefully. Failure to observe the instructions for use can endanger patient safety.

**Product Description**

SonoGuard is a sterile disposable tubular cover drape for ultrasound probes made of strong, thin polyethylene (LDPE) film that is telescopically folded and thus easily slipped onto the ultrasound probes to provide a protecting sheath over the transducer and its cable. To ensure the required ultrasound transmission from the sheathed transducer, coupling agents like gel are necessary inside the cover.

**Intended Use**

Sterile covers allow the use of non-sterile devices for aseptically performed diagnostic or therapeutic procedures on the human body. Correctly applied, the covers create a sterile barrier around the device. Equally, they prevent the contamination of probes by microorganisms or body fluids transmitting from the patient.

The covers are not intended for the use of the device in natural or surgical body orifices nor to prevent the device from being properly decontaminated between uses on patients.

**Contraindication and Potential Complications / Risks / Adverse Effects**

Contraindications are derived from the performed procedures and are to be taken from common scientific medical literature.

- A relative contraindication can be a known allergy against polyethylene.
- An accidental tear of the barrier film could lead to cross-contamination.

**Warnings and Precautions**

- For single patient, single procedure uses only. Do not re-sterilise and/or reuse it, as this can potentially result in compromised device performance and increased risk of inappropriate sterilisation and cross-contamination infection of the patient.
- Do not use covers, which have passed the expiration date.
- The package shall be checked to be intact before using. If the primary package is damaged, do not use the product.
- This product shall be used by qualified and experienced medical professionals only and any procedure must be conducted under strict sterile / aseptic conditions.
- Make sure that you have understood the intentions and functions of the product in order to use it safely and effectively.
- Unfold the cover softly; harsh pulling may in rare occasions lead to ruptures of the barrier film. Inspect visually for defects after unfolding.
- Reduce mechanical stress to the cover to a minimum. Accidental crushing, tearing or puncture of the cover by instruments, needles or equipment's accessories used may lead to failure of the barrier film.
- Do not remove or reposition adhesive strips and / or labels attached to the cover during use, as this may lead to failure of the barrier film.
- Do not bring the cover in contact with heated surfaces and heated instruments, which may lead to degradation of the barrier film.
- Dispose the cover after use in conformity with the institution's regulatory for contaminated devices and according to valid national regulations for such products.

**Application of SonoGuard (see Figures 1 - 6)**

1. The non-sterile person applies the ultrasound gel on the transducer face (Fig. 1).
2. The sterile person inserts the first hand into and the second hand over the protruding most inner pocket of the folded cover (Fig. 2).
3. The sterile person opens the second hand and pulls out the first hand, forming a tubular space in the most inner pocket for the transducer head (Fig. 3).
4. The sterile person keeps the second hand open to put the cover over the transducer presented by the non-sterile person (Fig. 4) and grabs it.
5. The sterile person extends the cover over the cable by pulling the turn-up (labelled) while holding the sheathed transducer firmly in the hand (Fig. 5).
6. Drape the cover tightly over the transducer face to remove any wrinkles and entrapped air. Fix the position of the cover by applying one of the adhesive strips or rubber bands supplied around the transducer (Fig. 6).
7. (Optional) With additional adhesive strip or rubber band provided, the loose-fit cover may be tightened behind the transducer to the cable to improve handling during use.

**Note**

- In case of pre-existing damage or wear of the transducer face, please do not use the cover.
- All serious incidents associated with the product should be reported to the manufacturer and the competent authority of the Member State where the user and / or the patient is established. The used cover must be preserved.

**Compatibility**

- The used materials of the cover are resistant to common disinfectants and coupling agents.
- In case a needle guide is applied to the sheathed transducer, make sure not to damage the film of the cover during use.

**eIFU**

The instructions for use can be downloaded at [www.promecon-medical.com](http://www.promecon-medical.com)

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PROMECON GmbH | Kreuzlerstraße 10 | 20095 Hamburg