
SonoGuard Plus

PROMECON

- DE** Sterile Abdeckung für Ultraschallsonden mit **ClearView** Haftfläche für den Schallkopf
- EN** Sterile cover for ultrasound probes with **ClearView** adhesive surface for the transducer
- FR** Housse stérile pour sondes à ultrason avec surface adhésive **ClearView** pour le transducteur
- IT** Copertura sterile per sonda ad ultrasuoni con superficie adesiva **ClearView** per il trasduttore
- ES** Funda estéril para sondas de ecografía con superficie adhesiva **ClearView** para el transductor
- PT** Capa de proteção esterilizada para sondas de ultrassom com superfície adesiva **ClearView** para o transdutor
- NL** Steriele hoes voor ultrasone probes met **ClearView**-kleefvlak voor de transducer
- DK** Steril afdækning til ultralydsprober med **ClearView** vedhæftningsflade til transduceren
- SE** Sterilt skydd för ultraljudssonder med **ClearView** självhäftande yta för transduktorn
- FI** Steriili suojus ultraäänisondeihin, jossa **ClearView**-tarrapinta ultraäänipäälle
- NO** Steril tildekning til ultralydsonder med **ClearView** hefteflate til lydhodet
- GR** Αποστειρωμένο κάλυμμα για κεφαλές υπερήχων με αυτοκόλλητη επιφάνεια **ClearView** για τον μετατροπέα

- DE** GEBRAUCHSANWEISUNG
- EN** INSTRUCTIONS FOR USE
- FR** MODE D'EMPLOI
- IT** ISTRUZIONI PER L'USO
- ES** INSTRUCCIONES DE UTILIZACIÓN
- PT** INSTRUÇÕES DE UTILIZAÇÃO

- NL** GEBUIKSAANWIJZING
- DK** BRUGSANVISNING
- SE** ANVÄNDARANVISNING
- FI** KÄYTTÖOHJE
- NO** BRUKSANVISNING
- GR** ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

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SonoGuard Plus

Sterile cover for ultrasound probes with ClearView adhesive surface for the transducer

Important

Please read the following information and instructions for use carefully. Failure to observe the instructions for use can endanger patient safety.

Product Description

SonoGuard Plus is a sterile disposable tubular cover drape for ultrasound probes made of strong, thin polyethylene (LDPE) film that is telescopically folded and thus easily slipped onto the ultrasound probes to provide a sheath over the transducer and its cable. On the inside, an adhesive surface, to be attached to the transducer face, ensures the required ultrasound transmission. There is no need to add coupling agents like gel inside the cover.

Intended Purpose

Sterile covers allow the use of non-sterile devices for aseptically performed diagnostic or therapeutic procedures on the human body. Correctly applied, the covers create a sterile barrier around the device. Equally, they prevent the contamination of the device by microorganisms or body fluids transmitting from the patient. The covers are not intended for the use of the device in natural or surgical body orifices nor to prevent the device from being properly decontaminated between uses on patients.

Contraindication and Potential Complications / Risks / Adverse Effects

- Contraindications are derived from the performed procedures and are to be taken from common scientific medical literature.
- A contraindication can be a known allergy against polyethylene.
- An accidental tear of the barrier film could lead to cross-contamination.

⚠ Warnings and Precautions

- For single patient, single procedure uses only. Do not re-sterilise and/or reuse it, as this can potentially result in compromised device performance and increased risk of inappropriate sterilisation and cross contamination infection of the patient.
- Do not use covers, which have passed the expiration date.
- The package shall be checked to be intact before using. If the primary package is damaged, do not use the product.
- This product shall be used by qualified and experienced medical professionals only and any procedure must be conducted under strict sterile / aseptic conditions.
- Make sure that you have understood the intentions and functions of the product in order to use it safely and effectively.
- Unfold the cover softly; harsh pulling may in rare occasions lead to ruptures of the barrier film. Inspect visually for defects after unfolding.
- Reduce mechanical stress to the cover to a minimum. Accidental crushing, tearing or puncture of the cover by instruments, needles or equipment's accessories used may lead to failure of the barrier film.
- Do not remove or reposition adhesive strips and / or labels attached to the cover during use, as this may lead to failure of the barrier film.
- Do not bring the cover in contact with heated surfaces and heated instruments, which may lead to degradation of the barrier film.
- Dispose the cover after use in conformity with the institution's regulatory for contaminated devices and according to valid national regulations for such products.

Application of SonoGuard Plus (see Figures 1 - 7)

1. The sterile person inserts the first hand into the protruding most inner pocket of the folded cover - the inside of the hand is facing the **ClearView** adhesive film - and spreads their fingers to prevent the formation of wrinkles on the adhesive surface (Fig. 1)
2. The sterile person removes the protective liner from the **ClearView** adhesive surface with the second hand. (Fig. 2)
3. The non-sterile person places the transducer on the **ClearView** adhesive surface, while the sterile person presses the adhesive surface against the transducer face from the inside. Then the first hand of the sterile person grabs the transducer head from the inside. (Fig. 3)
4. To sheath the transducer, the sterile person turns the folded cover inside out by pulling the flap (Label 1) with the second hand, while firmly holding the transducer with the first hand. (Fig. 4)
5. The non-sterile person extends the cover over the cable by pulling on the turn-up (Label 2) without touching the proximal end of the cover while the sterile person still holds the sheathed transducer firmly in the hand. (Fig. 5)
6. To remove any entrapped air or remaining wrinkles between the **ClearView** adhesive surface and the transducer face, the sterile person swipes over the related areas. (Fig. 6)
7. Fix the position of the cover by applying one of the adhesive strips or rubber bands provided around the transducer. (Fig. 7)
8. (Optional): With additional adhesive strip or rubber band provided, the loose-fit cover may be tightened behind the transducer to the cable to improve handling during use.

Note

- In case of pre-existing damage or wear of the transducer face, please do not use the cover.
- When removing the cover from the transducer, pull off the **ClearView** adhesive surface slowly and gently.
- All serious incidents associated with the product should be reported to the manufacturer and the competent authority of the Member State where the user and / or the patient is established. The used cover must be preserved.

Compatibility

- The used materials of the cover are resistant to common disinfectants and coupling agents.
- In case a needle guide is applied to the sheathed transducer, make sure not to damage the film of the cover during use.

eIFU

The instructions for use can be downloaded at: www.promecon-medical.com/ifu

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PROMECON GmbH | Kreuzlerstraße 10 | 20095 Hamburg