

Peters SURGICAL

Immeuble AURELIUM
1 Cours de l'Ile Seguin
92100 BOULOGNE-BILLANCOURT
FRANCE
Tél. +33 (0)1 48 10 62 62
Fax +33 (0)1 48 91 22 99
E-mail : peters@peters-surgical.com
www.peters-surgical.com

T'Lift® MD REF AW16280

A : Insertion / Insertion / Insertion / Inserimento/ inserción/ inserção/ indsættelse / innsetting / införing / invoeging

B : Retrait/ withdrawal/ Rückzug / ritiro / retirada / retirada / tilbagetrækning / uttak / tillbakadragande / opname

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INSTRUCTIONS FOR USE - En PETERS SURGICAL SINGLE-USE STERILE T'LIFT® RETRACTION SYSTEM

CAUTION

Please read these Instructions carefully prior to using the PETERS SURGICAL single-use T'Lift® retraction system. Failure to comply with the following recommendations may entail serious consequences. Laparoscopic surgery should only be performed by qualified doctors having been specifically trained for the purpose and having gained sufficient experience of use. All surgeons performing laparoscopic operations should be aware of the possible risks and complications involved in using such surgical technique; notably risks for the patient when manipulating instruments connected to an external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure. Patients, having a right to information according to Law, should be adequately briefed on the surgical procedure.

IMPORTANT!

- Warnings specific to a particular manipulation are stated throughout these Instructions (in bold type and preceded by ➢)
- These instructions solely contain information on the PETERS SURGICAL single-use T'Lift® and on how to use it. On no account should they be considered as an instructions' manual for using laparoscopic retraction systems in general.
- The PETERS SURGICAL T'Lift® is designed for single use only. On no account may it be used a second time over, whatever the reason or the means.
- The device is delivered in sterile packaging. The T'Lift® should never be resterilized, regardless of whether it has been used or not.
- Any reuse of this instrument shall affect its original sterility, may induce abnormal functioning or may lead to cross-contamination in patients. The manufacturer may not be held liable in any such a case. Reuse or resterilization may also compromise the integrity and functioning of the instrument, possibly causing injury, sickness, or even death of a patient.
- Responsibility with respect to the storage, the opening, the appropriate use, the handling, the insertion and the destruction of the PETERS SURGICAL T'Lift® lies solely with the operator. On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices for the storage, the opening, the appropriate use, the manipulation, the insertion and the destruction of the PETERS SURGICAL T'Lift®.
- Never use the sterile packaging and its contents in the following cases (whereupon everything should be destroyed):
 - Packaging severely crushed or clearly damaged;
 - Visible damp marks on the cardboard packaging or on the blister (shell or protective cover);
 - Perforation of the blister (shell or protective cover);
 - Expiry of the use-by date printed on the packaging.

INDICATIONS OF USE

The PETERS SURGICAL T'Lift® is an instrument exclusively designed for laparoscopic surgery. It should only be used by qualified practitioners having gained due experience thereof. The T'Lift® is designed to enable the suspension and presentation of anatomical sections during laparoscopic surgery.

CONTRAINDICATIONS

Use of the PETERS SURGICAL T'Lift® is strictly forbidden whenever a contraindication exists for laparoscopic surgery. On no account may PETERS SURGICAL be held liable for the assessment of any such contraindication for laparoscopic surgery. Not intended for use on new-born and infant.

DESCRIPTION

The T'Lift® is delivered in its packaging with:

- 2 T'Lift®,
- 2 locking clamps,
- 1 applicator (body and needle).

INSTRUCTIONS FOR USE

- 1- Check that the packaging (blister), serving to ensure sterile conditions for the medical instrument, is neither open nor damaged.
- 2- Check the use-by date printed on the packaging.
- 3- Open the sterile packaging according to proper surgical practices and set down the 5 elements within the operating zone.

- The PETERS SURGICAL T'Lift® and its packaging should be destroyed as of when sterility may no longer be guaranteed, notably in the following cases: expiry of the use-by date on the packaging; packaging severely crushed or clearly damaged; damp marks visible on the packaging; perforation of the blister (shell or protective cover).

4- Check that:

- the T'Lift® retraction system is in good condition;
- all components are duly present (T'Lift®, clamps and applicator).

- Any damaged T'Lift® should not be used.

5- Make sure, prior to surgery, that all laparoscopic instruments are adapted to the PETERS SURGICAL T'Lift®. Before inserting the PETERS SURGICAL T'Lift®, the operating zone, i.e. the abdominal wall, should be prepared as usual, in accordance with proper surgical practices.

- On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices when verifying compatibility of the laparoscopic instruments with the PETERS SURGICAL T'Lift®.

6- Remove the protection tube from the applicator; the user should be careful not to injure himself when handling the applicator (sharp object).

7- Insert the distal part of one of the T'Lift® into the main body of the applicator (see Plate A) and gently push until reaching the first marker.

8- Take hold of the T'Lift® fitted with the applicator and perforate the abdominal wall using proper surgical practices, while checking at the same time via endoscopy the needle's trajectory and exit point within the abdominal cavity.

9- Guided under endoscopy, grasp the anatomical section to be suspended, position it on the applicator's needle and pass right through the anatomical section.

10- Hold the body of the applicator and gently push the T'Lift® shaft until reaching the second marker.

11- Using the laparoscopic route, check the proper opening of the T'Lift® beneath the anatomical section chosen to be suspended.

12- Do not pull on the T'Lift® before having removed the applicator in order to avoid damaging it.

- Should the suspension operation fail (erroneous entry point, tissue rupture, unintentional opening of the instrument prior to tissue perforation, etc.), never reuse the T'Lift® in the case of it having been opened by way of the applicator, and remove it according to Plate B (diagram).

- If the T has been opened, do not attempt to withdraw it using the applicator.

13- Hold back the T'Lift® using an atraumatic clamp and remove the applicator via the external route, without releasing the anatomical section.

14- Place a clamp on the T'Lift® and block it in the desired position in order to correctly suspend the anatomical section. Do not "pull" the anatomical section using the T'Lift® only, but guide it along by way of an atraumatic clamp.

- Be sure to follow the regular surgical techniques, according to the surgical rules of the Art, in order to protect all vessels and organs from any possible injury.

15- Do not exert any undue tension on the anatomical section and check the quality of the tissues suspended.

- On no account may PETERS SURGICAL be held responsible in the case of failure by the operator to comply with the recommendations as regards surgical techniques. Residual risk remain, e.g. for vesical or vascular injuries, and tissue damage.

- Do not use any electrosurgical instruments directly in contact with the T'Lift® risk of damage to the T'Lift® and risk of burn .

16- When suspension is no longer required, follow the T'Lift® withdrawal procedure (see plate B): release the clamp of the T'Lift®; free up the distal part of the anatomical section of the T'Lift® via laparoscopic viewing; and take hold of the T'Lift® using a clamp (Figure ①). The T'Lift® should be held in place via its distal end (the largest and rounder part of the T'Lift®). Cut the T'Lift® on the outside of the abdominal cavity (Figure ②) and extract the distal part using the trocar via gentle pulling under laparoscopic control (Figure ③), thereby preventing damage to the tissues in the abdominal cavity or the T'Lift® being severed.

17- Immediately check both parts of the T'Lift® within the operation zone, in order to ensure that no parts are missing.

18- Repeat the above actions (from points 7 to 17) for the 2nd T'Lift® and the 2nd clamp.

Upon completion of surgery, the T'Lift® should be disposed of as operating-theatre waste:

- All T'Lift®, clamps and protective tubes should be placed in special containers reserved for soiled instruments,
- The applicator should be placed in the container reserved for sharp, traumatizing or perforating objects.

STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment adapted to sterile medical devices.
- For transportation, the above points should be observed.

Information for the United States (USA):

- US Federal Law restricts the sale of this instrument to doctors, or upon medical prescription only.
- Sterilized with ethylene oxide.
- Single use only - Not to be reused.

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