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En INSTRUCTIONS FOR USE OF DISPOSABLE STERILE TROCARS VECTEC

CAUTION

Read these instructions carefully before using VECTEC disposable trocars. Failure to respect the following recommendations may lead to serious consequences. Endoscopic operations must be carried out solely by doctors who have been specifically trained and have sufficient experience to use this technique. Any surgeon who undertakes an endoscopic operation must be aware of the risks and possible complications involved with this surgical technique (notably risks to the patient and to devices involving the use of any external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure). Ensure that the patient has been informed about the way the operation will take place, in accordance with current legislation.

IMPORTANT NOTE

- These instructions solely contain information concerning Vectec trocars and how to use them. They must in no way be considered as general instructions concerning techniques for the use of trocars.
- Vectec trocars are designed for single use only. They must in no way and on no account be used a second time.
- The trocars are delivered in sterile packaging. They must not be resterilised, regardless of whether or not they have been used.
- The reuse or resterilisation of this device compromises its sterility and may lead to infection or cross contamination between patients. Reuse or resterilisation can also compromise the integrity and functioning of the device and may lead to an injury or pathology or to the death of the patient.
- The warnings specific to a procedure have been grouped together with the present warning throughout this instruction sheet (in bold and preceded by the symbol: ).
- Responsibility for the storage, opening, suitable use, manipulation and disposal of the trocars lies solely with the operator. Vectec shall in no case be held responsible if the rules for good surgical practice are not applied by the operator with respect to the storage, opening, suitable use, manipulation and disposal of the trocar. The manufacturer is likewise relieved of all responsibility in the case of reuse and/or resterilisation.

FIELD OF USE

Vectec trocars are instruments designed exclusively for endoscopic surgery. They allow the passage of endoscopic instruments.

COUNTER INDICATIONS

The use of Vectec trocars is strictly excluded when there is a counter indication for endoscopic surgery.

- Vectec can in no case be held responsible for the assessment of any counter indication for endoscopic surgery.

DESCRIPTION OF THE TROCAR

Vectec trocars are comprised of the trocar chuck (part ② figure A and part ③ figure B) and the threaded radio-transparent trocar sleeve (part ① figures A and B).

The inside diameter of the sleeve is 5.5 mm for references 5 053/5 053H/5 053T/5 068/5 100, 12 mm for references 12 060/12 100/12 100H/11 100S/11 060/11 100.

The length of the trocar varies according to the reference and is marked on the product packaging.

A reducing device is supplied for references 12 060/12 100/12 100H/11 100S/11 060/11 100 (part ④ in figure B) for use with instruments of 5 mm in diameter.

The trocar sleeve is fitted with a flexible valve, to keep gas leaks to a minimum.

There is a compatible Luer-Lock connection for insufflation:

- With a stopper/roref. 5 068/5 100 (see figure A),
- With a Luer-Lock valve for ref. 12 060/12 100/12 100H/11 100S/11 060/11 100 (see figure B).

UTILISATION

Preparation

1. Open the multipack box and take out the medical device,
2. Check that the packaging (blister) ensuring sterile conditions for the medical device is neither open nor damaged.
3. Check the expiry date for use marked on the package.
4. Open the sterile packaging according to proper professional practice.

➤ Vectec trocars and their packaging must be destroyed if sterility is not guaranteed; notably in the following cases: if the date marked on the packaging has been exceeded, the packaging has been severely crushed or is clearly damaged, traces of humidity are visible on the packaging, the blister is punctured (bubble or membrane seal).

5. Check:

- the condition and proper functioning of the trocar, sleeve and chuck,
- the presence of components: sleeve, chuck and, depending on the references, the reducing device and the fixing cone (see figures A and B).

➤ Damaged trocars must not be used.

6. Before starting a procedure, check that all the endoscopic instruments are suitable for use with Vectec trocars (diameters and lengths; see § "Description"). Check the condition and proper functioning of the instruments before inserting them in the trocars, particularly that of the insulation sheath.

➤ To limit leakage and avoid loss of pneumoperitoneum, do not use instruments with a diameter of less than 4.2 mm with trocars of 5 mm or of less than 8.5 mm with trocars of 12 mm (4.2 with the 5 mm reducing device), see figures F and G.

➤ The use of instruments with maximal diameter can limit the flow of insufflation. In this case, it is recommended to insufflate through another trocar.

➤ In no case can Vectec be held responsible if the recommendations of good surgical practice have not been applied by the operator to check the compatibility (diameters and lengths) of the endoscopic instruments with Vectec trocars.

7. Before inserting the Vectec trocar, the operating field – the abdominal wall – must be prepared according to the rules of good surgical practice. Before inserting the trocar, verify the anatomical disposition of the

patient in order to minimise any risk of lesion, either vascular or of the internal organs.

Inserting the trocar

8. Insert the chuck at the centre of the conical opening of the trocar sleeve.

➤ Insert the cutting tips carefully to avoid harming the joints and causing a leak.

9. Testing the shield chuck (trocar ref. 11100S): Pick up the trocar, with the chuck and sleeve assembled. Check that the shield moves freely (see figure D: mobility up and down the trocar) if not, refer to § "Resetting the shield". The trocar is now ready for insertion in the abdominal wall.

10. Resetting the shield (trocar ref. 11100S): Before inserting the trocar into the abdominal wall, the shield must be set after inserting the chuck in the sleeve. At this stage, the shield is in the locked position, so with the safety catch on. The shield covers the 3 mm tip and can no longer move back up the triggering groove. To reset the shield, take the assembled chuck and sleeve in one hand, with the tip pointing outwards. Grasp the shield with the other hand and turn the shield anticlockwise. Rotate through 180° i.e. half a turn, until it clicks into the set working position (see figure E).

➤ The locked shield stops the tip from damaging the seal around the trocar sleeve. Setting and testing should be done gently to avoid any injury from the tip. Setting and testing do not need any great effort or "forcing", and the shield has only one possible operational movement. Take care: the chuck is single use only so the cutting end becomes blunt with use. The considerable stresses placed on the shield only allow it to be triggered 3 or 4 times.

11. Place the trocar in a sterile area.

➤ Do not damage the tip of the chuck due to a shock or similar movement.

12. Select the point of insertion. Make the incision according to the usual technique used in endoscopic procedures in compliance with the rules of surgical practice. Trocars ref. 12 100H must be used according to "Open laparoscopy" techniques [secure the position of the trocar in the abdominal wall using a fixing cone (part ④ in figure B) and a purse-string suture].

13. Insert the trocar with a slight alternating, rotating motion, under endoscopic control, exerting constant pressure on the abdominal wall (for versions 5 053/5 053H/5 068/5 100/12 060/12 100/11 060/11 100).

➤ Attention: trocar 11 100S is inserted WITHOUT rotation.

➤ Apply the habitual techniques according to the rules of surgical practice for protecting the vessels and organs from possible injuries, particularly in the case of the first trocar which cannot be inserted under endoscopic control.

14. For the first trocar, remove the chuck and use the endoscope to check the integrity of the abdominal wall and the correct placement of the sleeve.

15. Screw the trocar in clockwise until the trocar sleeve is fixed perfectly straight and safely in the abdominal wall. For trocars inserted under endoscopic control, remove the chuck from the sleeve with a regular motion.

Insertion of instruments

16. Insert and remove endoscopic instruments of a suitable diameter (see § «Description») through the trocar sleeve.

➤ Insert and remove instruments carefully, particularly instruments with a cutting tip, to avoid any damage to sealing joints which may lead to pneumoperitoneum loss.

End of procedure

17. Smoothly unscrew the sleeve from the trocar in an anticlockwise direction, to extract it from the abdominal wall. Perform suture according to the habitual techniques of surgical practice. After the procedure, the trocar must be handled as operating theatre waste and the chuck must be disposed of in a special container for cutting or perforating objects that might cause injury.

WARNING

Do not use but destroy the sterile package together with its contents in the following situations:

- Packaging is severely crushed or clearly damaged.
- Traces of humidity are visible on the packaging.

STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment that is suitable for sterile medical devices.
- For transport conditions apply the above points.

Information for the United States (USA):

- American federal law limits the sale of this device to doctors or with a doctor's prescription.
- Sterilized by ethylene oxide.
- Single use only – do not reuse.

In the case of any adverse event (or materiovigilance), kindly contact Vectec at this address: regulatory@vectec.fr

First CE mark: MARCH 1999.

Manufactured by:

VECTEC

Bioparc - 03270 HAUTERIVE - FRANCE

Tel. : + 33 (0)4 70 58 41 70

Fax : + 33 (0)4 70 58 41 79

vectec@vectec.fr

REF

VEC-1110S

VEC-12100 H

VEC-5053

VEC-5068

VEC-5100

VEC-12060

VEC-12100

CE 0459

STERILE EO

DEHP