

Cover for surgical light handle

Sterile cover for surgical lights

- o Sterile cover for surgical light handle
- o Minimizes contamination risk of surgical field
- o Universal for all handles*
- o Latex free, DHP free, tack-free
- o Single or double sterile packacking

* Handles with 20-40 mm diameter
Outer diameter of fixing disc 11.8cm



Order code	Dimensions	SU
OP.00.LHC-01	Pouch: 9.5x13cm (for handles with Ø 20-40mm)	100 / box (single packaging)
OP.00.LHC-02	Pouch: 9.5x13cm (for handles with Ø 20-40mm)	50 x 2 / box (double packaging)
OP.00.LHC-03	Pouch: 9.5x15cm (for handles with Ø 20-40mm)	100 / box (single packaging)

